



EMERALD HEALTHCARE

WE CARE, SO YOU CAN TOO

HOW TO PROGRAM OUR NEW BED OR CHAIR PAD

NOTE: Before using, you will need to program your new pads to your alarm. Taking a small paperclip, find the on/off button on the transmitter on the new pad and push down on this button, you should see the red status light blink a couple of times then turn off.

1. Have your pad laid out in front of you and have your alarm ready in the chime position.
2. On the back of the ALARM monitor is a hole with the word “Program” below it. Taking a very slim screwdriver or paperclip, insert into that hole, you will feel a button inside. Push down on the button and hold down – the two outer lights on the alarm monitor will flash and then stay on. Remove the screwdriver/paperclip. (if the lights blink and go off, try again, they need to stay on once the pen/screwdriver is removed)
3. Now apply pressure to the pad for about 5 seconds, after this time release the pressure.
4. Quickly taking the screwdriver/paperclip, push into the program hole on the back of the alarm again. Push down on the button inside and then release. The lights on the front of the alarm will flash and then go off, your alarm should also chime.

To Test: After programming put your alarm into the off mode using the slide switch on the side, wait a few seconds for the alarm to reset and then place it back into chime mode. Press down on the pad for a few seconds and then release, your alarm should chime. If you do not get a chime you need to go through this process again.

If you require any further information or assistance, please do not hesitate in contacting us.



EMERALD HEALTHCARE
WE CARE, SO YOU CAN TOO

T: 01 4245103

E: info@emeraldhealthcare.ie