Prevent Falls & Wandering



Safety Auto-Reset

TL-2100S

Tamperproof Monitor Can Only be Turned Off by Caregiver





Automatically Turns On When Resident Sits or Lies on Pad

Simply plug a sensor pad into the alert and place the weight-sensing pad under the resident in a chair or bed. When the resident gets up and weight is removed, the alert will sound notifying caregivers that the resident is on the move. This alert features an Auto-Reset® function which will automatically reset and silence the alert when weight is re-applied to the sensor pad.

Optional Settings

When set to TamperProof®, the alert can ONLY be reset using the provided Reset Key. Selecting this setting ensures that only the designated caregiver can silence the alert.

Key Features:

- Full-Range Volume Contro
- Nurse Call Capability enables use with existing hard wired purse call systems.
- Optional Setting: TamperProof® Reset Ke allows control of who is able to reset the
- Batteries not included (3 x AA
- Optional AC adapter (AC-05)
- Protective Boot

For use with one of the following Corded Sensors:

PPB-WI - Bed Sensor Pad 20" x 30"

PPB-RI — Bed Sensor Pad 10" x 30"

PPC-WI – Chair Sensor Pad 10" x 15"

TL-2109 - Easy Release Seat Belt 54"

TL-2109V — Quick Release Seat Belt w/hook & loop 44"

FM-05 – Floor Mat 24" x 36" (gray)

FM-06 – Floor Mat 20" x 30" (under rug) (gray)

FM-07 – Floor Mat 24" x 48" (gray)

WC-02 — Weight-Sensing Fall Mattress 29.5"x 70.5"x 1.5" (blue)

LM-01 — Weight-Sensing Impact Landing Mat 24"x71"x1" (gray)

Tel: 01 - 424 5103 Email: info@emeraldhealthcare.ie



433-EC

Low Cost Fall Prevention





No Cords! Place monitor where you need it.

Low Cost Solution

The Economy CordLess Fall Prevention Alert provides affordable and effective cord free monitoring with a range up to 300 ft. Ideal for home or small facility use.

Easy-To-Use

The Economy CordLess Fall Prevention Alert is easy-to-use and works with all CordLess Sensor Pads, Floor Mats, Nurse Call Buttons, and Motion Sensors. Monitor is light-weight and can be programmed with up to six (6) wireless components to serve as a portable caregiver alert.

Remove Noise From the Room

Alert can be mounted away from the bed, removing in-room alert noise. Adjustable volume and gentle chime alert won't startle residents and offers QUIETER fall prevention.

- Compliant for Non-Restraint
- Nurse Call Capability enables use with

#1 Smart Choice

For use with the following CordLess Sensors: See Images Next Page

433-NC - Nurse Call Button

433-TB - Touch Call Button

433-RB - Remote Reset Button

433-MS - Motion Sensor

433-EXT - Window/Door Exit Sensor

GBT-WI - CordLess Bed Pad 20" x 30"

GBT-RI - CordLess Bed Pad 10" x 30"

GCT-WI - CordLess Chair Pad 10" x 15"

FMT-05C - CordLess Floor Mat 24" x 36" (grey)

FMT-06C - CordLess Floor Mat 20" x 30" (grey) (under rug)

FMT-07C - CordLess Floor Mat 24" x 48" (grey) LM-01C - Cordless Landing Mat 24" x 72"

Tel: 01 - 424 5103

Email: info@emeraldhealthcare.ie

